

# scRGOT Coder Upgrade: Spatial Transcriptomics




Corinne Strawser, PhD


Bioinformatics Scientist

Dr. Elaine Mardis & Dr. Katie Miller Labs

Institute for Genomic Medicine



# Why spatial?



Complex Tissue

*Networks of many cells each with individual gene expression patterns that regulate interactions and ultimately tissue function*



Bulk RNA-seq

*tissue-average gene expression*



Single-cell RNA-seq


*Cell-specific gene expression*




Spatial Transcriptomics

*region-, cell-, or subcellular-specific gene expression with spatial context*

# Spatial transcriptomics has broad applications





# Spatial transcriptomics is not “new”...







# ...but it's expanding!





# 10X Genomics Visium


- Visium slides has two or four capture areas of either 6.5 x 6.5 mm or 11 x 11 mm.
- ~5,000 capture spots (6.5 x 6.5 mm capture area)
- ~14,000 barcoded spots (11 x 11 mm capture area)
- Each spot has millions of oligonucleotides with:
  - poly(dT) sequence for the capture of polyadenylated molecules
  - UMI used for the identification of duplicate molecules that arise during the library preparation and sequencing process
  - Spatial Barcode, which is shared by all oligonucleotides within each individual gene expression capture spot






# 10X Genomics Visium HD



- The HD slide has two capture areas, each 6.5 x 6.5 mm
- ~ 11 million 2 x 2  $\mu\text{m}$  barcoded squares without gaps
- Oligonucleotides each have
  - poly(dT) sequence for capture of the ligation product
  - UMI for the identification of duplicate molecules that arise during the library preparation and sequencing process
  - Spatial Barcode



Visium HD Gene Expression Probe-Based Library



# 10X Genomics Visium Workflows



# 10X Genomics Visium and Visium HD


## PRO

- Very well documented
- Microscopy image
- FFPE, fresh frozen, and fixed frozen (HD and LD v2) or fresh frozen (LD v1)
- Protein co-detection
- Customizable
- Any species (LD v1)

## CON

- U-Pick: subcellular (HD) or multi-cellular (LD)
- Probe-based (HD & LD v2) – *Visium HD 3' coming soon*
- CytAssist required (HD & LD v2)
- Mouse or human (HD & LD v2)

# Spatial Transcriptomics Analysis



# spaceranger mkfastq

## mkfastq sample sheet

```
$ run_spaceranger_mkfastq.sh samplesheet.csv ×
igm > projects > 241014_GSL-LPRV-4138_Visium > fastq > samplesheet.csv > data
1 Lane, Sample, Index,
2 *,NP-20-0025854_C1_Pre_Gene,SI-TS-E6,
3 *,NP-20-0043174_D4_Post_Gene,SI-TS-A5,
4 *,NP-20-0025854_C1_Pre_Protein,SI-NT-E1,
5 *,NP-20-0043174_D4_Post_Protein,SI-NT-F1,
6 *,NP-19-0035324_A1_Gene,SI-TS-D7,
7 *,NP-19-0042829_B3_Gene,SI-TS-D8,
8 *,NP-19-0035324_A1_Protein,SI-NT-C1,
9 *,NP-19-0042829_B3_Protein,SI-NT-D1,
```

## mkfastq script

```
$ run_spaceranger_mkfastq.sh × samplesheet.csv
igm > projects > 241014_GSL-LPRV-4138_Visium > fastq > $ run_spaceranger_mkfastq.sh
1 #SBATCH -q all.q
2 #SBATCH -j y
3 #SBATCH -cwd
4 #SBATCH -S /bin/bash
5 #SBATCH -terse
6 #SBATCH -pe smp 16
7 #SBATCH -N cellranger_mkfastq
8
9 module load spaceranger_3.0.1
10
11 spaceranger mkfastq \
12 --samplesheet=samplesheet.csv \
13 --run=/igm/runs/IGM_Seq10/241018_A00498_0918_BHH27NDRX5 >& bclconv.out
```

\* note you'll have to change the header if you're using Slurm as your job scheduler (I'm using SGE here).

# spaceranger count samplesheet

```
NP-19-0035324_A1.csv × $ NP-19-0035324_A1.sh ●  
igm > projects > 241014_GSL-LPRV-4138_Visium > counts > NP-19-0035324_A1 > NP-19-0035324_A1.csv > data  
1 fastqs,sample,library_type  
2 /igm/projects/241014_GSL-LPRV-4138_Visium/fastq/HH27NDRX5/outs/fastq_path/HH27NDRX5,NP-19-0035324_A1_Gene,Gene Expression  
3 /igm/projects/241014_GSL-LPRV-4138_Visium/fastq/HH27NDRX5/outs/fastq_path/HH27NDRX5,NP-19-0035324_A1_Protein,Antibody Capture
```




# spaceranger count script

```
NP-19-0035324_A1.csv  NP-19-0035324_A1.sh •
igm > projects > 241014_GSL-LPRV-4138_Visium > counts > NP-19-0035324_A1 > $ NP-19-0035324_A1.sh
1  ## -j y
2  ## -cwd
3  ## -S /bin/bash
4  ## -N spaceranger_count
5  ## -q all.q
6  ## -pe smp 16
7  ## -V
8
9  # set up
10 dir="/igm/projects/241014_GSL-LPRV-4138_Visium/counts/NP-19-0035324_A1"
11 id="NP-19-0035324_A1_Pre"
12 libraries="/igm/projects/241014_GSL-LPRV-4138_Visium/counts/NP-19-0035324_A1/NP-19-0035324_A1.csv"
13 image="/igm/projects/240725_GSL-LPRV-4007/images/H&E/NP-19-0035324_A1_Pre_transformed.jpeg"
14 cytaimage="/igm/projects/240725_GSL-LPRV-4007/images/assay_CAVG10321_2024-07-18_11-19-22_V53A24-007_1721316809_CytAssist/
CAVG10321_2024-07-18_11-33-29_2024-07-18_11-19-22_V53A24-007_A_NP-19-0035324_A1.tif"
15 alignment="/igm/projects/240725_GSL-LPRV-4007/images/assay_CAVG10321_2024-07-18_11-19-22_V53A24-007_1721316809_CytAssist/
final_alignment_np_19_0035324_a1.json"
16 transcriptome="/igm/home/cnh008/ref/C134_M032/refdata-gex-GRCh38-2020-A-C134"
17 probes="/igm/home/cnh008/ref/C134_M032/Visium_Human_Transcriptome_Probe_Set_v2.0_GRCh38-2020_M032_C134_IDT_custom_probes.csv"
18 protein="/igm/apps/10X_chromium/spaceranger-3.0.0/feature_refs/Visium_Human_Immune_Cell_Profiling_Panel_v1.0.csv"
19
20 # script
21 cd ${dir}
22 module purge
23 module load spaceranger_3.0.1
24
25 spaceranger count \
26   --id=${id} \
27   --transcriptome=${transcriptome} \
28   --libraries=${libraries} \
29   --probe-set=${probes} \
30   --image=${image} \
31   --cytaimage=${cytaimage} \
32   --create-bam=true \
33   --loupe-alignment=${alignment} \
34   --feature-ref=${protein} \
35   --localcores=16 \
36   --localmem=64
```

\* note you'll have to change the header if you're using Slurm as your job scheduler (I'm using SGE here).

# Raw data processing:



# Downstream analysis: choose your own adventure!



:



Giotto



squidpy

ST  
LEARN



:




Giotto



(formerly STUtility)



# Today's data: Spatial transcriptomic map of aging mouse brain




# Regions of Interest (ROI)

## Mouse, Adult, 3D Coronal

Acronym: Isocortex  
Name: Isocortex


Adult, 3D Coronal Atlas



## Mouse, Adult, 3D Coronal

Acronym: HIP  
Name: Hippocampal region

Adult, 3D Coronal Atlas



# Other spatial -omics

- **DNA sequencing:** slide-DNA-seq (*Zhao et al., 2022*)
- **Epigenomics:** ATAC-seq (*Deng et al., 2022*) and CUT&TAG (*Deng et al., 2022*) both available with FlowGel assay by AtlasXomics
- **Proteomics:** CODEX (*Goltsev et al., 2018*), Immuno-SABER (*Saka et al., 2019*), scDVP (*Rosenberger et al., 2023*), RIBOmap (*Zeng et al., 2023*), CITE-seq (*Liu et al., 2023*), Stereo-CITE-seq (*Liao et al., 2023*)
- **Metabolomics:** MSI (*Unsihuay et al., 2021*) and scSpaMet (*Hu et al., 2023*)
- **Temporally-resolved:** intMEMOIR (*Chow et al., 2021*), Spatial iTracer (*He et al., 2022*), Space-TREX (*Ratz et al., 2022*)
  - Enable lineage tracing using CRISPR-Cas9 scarring or lentiviral barcoding!



# Experimental Considerations

## 1. Experimental design

- Hypothesis generation or testing?

## 2. Target (gene or cell) detectability


- Rare cells will be hard to identify with a low-res technology
- Transcripts with low expression may be difficult to profile with a technology that has low capture efficiency

## 3. Tissue type and species

## 4. Availability of reference data

- For cell deconvolution

## 5. Tissue section type and quality



# What challenges remain?

- Compatibility with wider range of tissues
- Optimize data storage and processing speed
- Optimize spatial resolution and capture efficiency
- Reduce signal diffusion
- ST-specific normalization and batch correction
- Registration across slices
- Integration with multi-omics
- Increased throughput, reduced cost, greater accessibility

